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
## Bibliografia

- <sup>1</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1993). Periodontal regeneration of human infrabony defects. I. Clinical Measures. *Journal of Periodontology* 64, 254-260.
- <sup>2</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1993). Periodontal regeneration of human infrabony defects. II. Re-entry procedures and bone measures. *Journal of Periodontology* 64, 261268.
- <sup>3</sup> Tonetti, M., Pini-Prato, G. & Cortellini, P. (1993). Periodontal regeneration of human infrabony defects. IV. Determinants of the healing response. *Journal of Periodontology* 64, 934-940.
- <sup>4</sup> Tonetti, M., Pini-Prato, G. & Cortellini, P. (1995). Effect of cigarette smoking on periodontal healing following GTR in infrabony defects. A preliminary retrospective study. *Journal of Clinical Periodontology* 22, 229-234.
- <sup>5</sup> Tonetti, M., Pini-Prato, G. & Cortellini, P. (1996). Factors affecting the healing response of intrabony defects following guided tissue regeneration and access flap surgery. *Journal of Clinical Periodontology* 23, 548-556.
- <sup>6</sup> Cortellini, P. & Tonetti, M.S. (2000). Focus on intrabony defects: guided tissue regeneration (GTR). *Periodontology* 2000 22, 104-132.
- <sup>7</sup> Cortellini P, Tonetti MS. (2005). Clinical performance of a regenerative strategy for intrabony defects: scientific evidence and clinical experience. *J Periodontol.* Mar;76(3):341-50.
- <sup>8</sup> Cortellini P. Reconstructive periodontal surgery: a challenge for modern periodontology. *Int Dent J* 2006; 56 suppl 1:250-255.
- <sup>9</sup> Cortellini, P., Carnevale, G., Sanz, M. & Tonetti, M.S. (1998). Treatment of deep and shallow intrabony defects. A multicenter randomized controlled clinical trial. *Journal of Clinical Periodontology* 25, 981-987.
- <sup>10</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). The modified papilla preservation technique with bioresorbable barrier membranes in the treatment of intrabony defects. Case reports. *International Journal of Periodontics and Restorative Dentistry* 14, 8-15.
- <sup>11</sup> Cortellini, P., Prato, G.P. & Tonetti, M.S. (1999). The simplified papilla preservation flap. A novel surgical approach for the management of soft tissues in regenerative procedures. *International Journal of Periodontics and Restorative Dentistry* 19, 589-599.
- <sup>12</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1995). Periodontal regeneration of human infrabony defects with titanium reinforced membranes. A controlled clinical trial. *Journal of Periodontology* 66, 797-803.
- <sup>13</sup> Cortellini, P. & Tonetti, M. (1999). Radiographic defect angle influences the outcome of GTR therapy in intrabony defects. *Journal of Dental Research* 78, 381 (abstract).
- <sup>14</sup> Tonetti MS, Cortellini P, Lang NP, Suvan JE, Adriaens P, Dubravec D, Fonzar A, Fourmouis I, Rasperini G, Rossi R, Silvestri M, Topoll H, Walkkamm B, Zybutz M. (2004). Clinical outcomes following treatment of human intrabony defects with GTR/bone replacement material or access flap alone. A multicenter randomized controlled clinical trial. *J Clin Periodontol.* Sep;31(9):770-6.
- <sup>15</sup> Cortellini, P., Tonetti, M.S., Lang, N.P., Suvan, J.E., Zucchelli, G., Vangsted, T., Silvestri, M., Rossi, R., McClain, P., Fonzar, A., Dubravec, D. & Adriaens, P. (2001). The simplified papilla preservation flap in the regenerative treatment of deep intrabony defects: clinical outcomes and postoperative morbidity. *Journal of Periodontology* 72, 1701-1712.
- <sup>16</sup> Linares A, Cortellini P, Lang NP, Suvan J, Tonetti MS. Guided tissue regeneration/deproteinized bovine bone mineral or papilla preservation flap alone for treatment of intrabony defects. II: radiographic predictors and outcomes. *J Clin Periodontol* 2006; 33: 351-358.
- <sup>17</sup> Tonetti MS, Fourmouis I, Suvan J, Cortellini P, Bragger U, Lang NP; European Research Group on Periodontology (ERGOPERIO). (2004). Healing, post-operative morbidity and patient perception of outcomes following regenerative therapy of deep intrabony defects. *J Clin Periodontol.* Dec;31(12):1092-8.
- <sup>18</sup> Cortellini, P. & Tonetti, M.S. (2001). Microsurgical approach to periodontal regeneration. Initial evaluation in a case cohort. *Journal of Periodontology* 72, 559-569.
- <sup>19</sup> Mayfield L, Tonetti MS, Cortellini P, Lang NP. Microbial colonisation pattern predict the outcomes of surgical treatment of intrabony defects. *J Clin Periodontol* 2006; 33:62-68.
- <sup>20</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1994). Periodontal regeneration of human infrabony defects. V. Effect of oral hygiene on long term stability. *Journal of Clinical Periodontology* 21, 606-610.
- <sup>21</sup> Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). Long term stability of clinical attachment following guided tissue regeneration and conventional therapy. *Journal of Clinical Periodontology* 23, 106-111.
- <sup>22</sup> Cortellini P, Tonetti MS. (2004). Long-term tooth survival following regenerative treatment of intrabony defects. *J Periodontol.* May;75(5):672-8.

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- Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). Periodontal regeneration of human intrabony defects with bioresorbable membranes. A controlled clinical trial. *Journal of Periodontology* 67, 217-223.
- Cortellini, P., Pini-Prato, G. & Tonetti, M. (1995). The modified papilla preservation technique. A new surgical approach for interproximal regenerative procedures. *Journal of Periodontology* 66, 261-266.
- Cortellini, P. & Tonetti, M. (2000). Evaluation of the effect of tooth vitality on regenerative outcomes in intrabony defects. *Journal of Clinical Periodontology* 28, 672-679.
- Cortellini P, Tonetti MS (2007) A minimally invasive surgical technique (MIST) with enamel matrix derivate in the regenerative treatment of intrabony defects: a novel approach to limit morbidity. *Journal of Clinical Periodontology*; 34: 87-93.
- Cortellini P, Tonetti MS (2007) Minimally invasive surgical technique (MIST) and enamel matrix derivative (EMD) in intrabony defects. (I) Clinical outcomes and morbidity. *Journal of Clinical Periodontology*; 34: 1082-1088.
- Tonetti, M. S., Pini-Prato, G. P., Williams, R. C. & Cortellini, P. (1993). Periodontal regeneration of human infrabony defects. III. Diagnostic strategies to detect bone gain. *Journal of Periodontology* 64, 269-277.
- Tonetti, M., Pini-Prato, G. & Cortellini, P. (1996). Guided tissue regeneration of deep intrabony defects in strategically important prosthetic abutments. *International Journal of Periodontics and Restorative Dentistry* 16, 378-387.
- Tonetti, M., Cortellini, P., Suvan, J.E., Adriaens, P., Baldi, C., Dubravec, D., Fonzar, A., Fourmouis, I., Magnani, C., Muller-Campanile, V., Patroni, S., Sanz, M., Vangsted, T., Zabalegui, I., Pini Prato, G. & Lang, N.P. (1998). Generalizability of the added benefits of guided tissue regeneration in the treatment of deep intrabony defects. Evaluation in a multi-center randomized controlled clinical trial. *Journal of Periodontology* 69, 1183-1192.
- Tonetti, M., Lang, N.P., Cortellini, P. et al. (2002). Enamel matrix proteins in the regenerative therapy of deep intrabony defects. A multicenter randomized controlled clinical trial. *Journal of Clinical Periodontology*;29:317-25.
- Tsitoura E, Tucker R, Suvan J, Laurell L, Cortellini P, Tonetti M. Baseline radiographic defect angle of the intrabony defect as a prognostic indicator in regenerative periodontal surgery with enamel matrix derivative. *J Clin Periodontol* 2004; 31:643-647.

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©Geistlich Pharma AG